PHOENIX ELECTRODE COMPANY CYANIDE ION ELECTRODES INSTRUCTION MANUAL

GENERAL INSTRUCTIONS

Introduction

The pHoenix Cyanide Ion Electrodes are used to quickly, simply, acc urately, and economically measure cyanide ions in aqueous solutions.

IMPORTANT NOTICE

Acidic cyanide solutions produce hydrogen cyanide (HCN) gas, highly toxic whether breathed or absorbed through the skin. Use of the p roper and recommended ionic strength adjuster (ISA) will keep the s olution pH above 10. If measurements in acidic solution are necess ary (decomplexing procedure as given in the section **Complexation**), **THE PROCEDURE MUST BE DONE IN A HOOD**.

Use a pipet bulb when pipeting cyanide solutions, as these solution s are highly toxic.

Required Equipment

- 1. ApH/mV meter or an ion meter, either line operated or portable.
- 2.

Semi-logarithmic 4-cycle graph paper for preparing c alibration curves when using the meter in the mV mode.

3. A magnetic stirrer.

4.

The pHoenix Cyanide Ion Electrode, Cat. No. CN01501 (reference electrode necessary), the pHoenix Cyanide Ion Combin ation Glass Electrode, Cat. No. CN01502, or the pHoenix Cyanid e Ion Combination Epoxy Electrode, Cat. No. CN01503.

5.

The pHoenix Double Junction Reference Electrode, Cat. No. 5731429 (for use with the CN01501) with pHoenix Filling So lution, Cat. No. R001011, in the inner junction and with pHoen ix Filling Solution, Cat. No. R001015, in the outer junction.

6. Lab-ware made of plastic, not glass.

7.

pHoenix Polishing Paper, Cat. No. 948201, to polish dirt y or etched electrode membranes.

Required Solutions

1.

Deionized or distilled water for solution and standard preparation.

2.

pHoenix Cyanide Ionic Strength Adjuster (ISA), 10M NaOH, Cat. No. CNOISO1. To prepare the ISA from your own laborator y stock, fill a 1000 ml beaker with about 900 ml of distilled water. While gently stirring the solution under a hood, slowl y add 400 grams of reagent-grade sodium hydroxide. Transfer t he solution quantitatively to a one liter volumetric flask aft er the solid NaOH has dissolved and the beaker has cooled. Di lute to the mark with distilled water, cap, and upend several times to thoroughly mix the solution. Store in a plastic bottle.

3.

Cyanide Standard, $1X10^{-2}M$. To prepare this solution from your own laboratory stock, add 10 ml of ISA and about 500 ml o f distilled water to a one liter volumetric flask. Add 0.49 g rams of dry, reagent-grade sodium cyanide, NaCN, to the soluti on and swirl the flask gently to dissolve the solid. Dilute t o the mark with distilled water, cap, and upend the flask seve ral times to thoroughly mix the contents. Store all standards in plastic bottles and prepare weekly.

4.

Cyanide Standard, 1000 ppm. To prepare this solution fro m your own laboratory stock, add 10 ml of ISA and about 500 ml of distilled water to a one liter volumetric flask. Add 1.88 grams of dry, reagent-grade NaCN and swirl the flask gently t o dissolve the solid. Dilute to the mark with distilled water , cap, and upend the flask several times to thoroughly mix the contents. Store all standards in plastic bottles and prepare weekly.

GENERAL PREPARATION

Electrode Preparation

Remove the rubber caps covering the electrode tips and the rubber i nsert covering the filling hole of the cyanide combination ion elec trode or the reference electrode. Fill the reference electrode or the combination electrode with the filling solution shipped with th e electrode to a level just below the fill hole. No preparation is required with a sealed reference electrode. Connect the electrode s to the proper terminals as recommended by the meter manufacturer.

Electrode Slope Check (with pH/mV meter) (check electrodes each day)

1.

To a 150 ml beaker, add 100 ml of distilled water and 1 ml of 10M NaOH ISA. Place the beaker on a magnetic stirrer an d begin stirring at a constant rate. After assuring that the meter is in the millivolt mode, lower the electrode tips into the solution.

2.

Using a pipet, add 1 ml of $1.0 \times 10^{-2} M$ or 1000 ppm standard to the beaker. When the reading is stable, record the mV reading.

3.

Using a pipet, add 10 ml of the same standard used above to the beaker. After the reading has stabilized, record the mV reading.

4.

The electrode is operating correctly if the mV potential has changed by 57+2 mV, assuming the solution temperature is between 20° and $25^{\circ}C$. See the **TROUBLESHOOTING** sections if the potential change is not within this range.

<u>Slope</u> is defined as the change in potential observed when the conce ntration changes by a factor of 10.

Electrode Slope Check (with an ion meter) (check electrodes each day)

1.

Prepare standard cyanide solutions whose concentrations vary by tenfold. Use either the $1.0 \times 10^{-2} M$ or the 1000 ppm cyan ide standard. Use the serial dilution method for this preparation.

2.

To a 150 ml beaker, add 100 ml of the lower value standa rd and 1 ml of ISA. Place the beaker on a magnetic stirrer an d begin stirring at a constant rate. Lower the electrode tips into the solution. Assure that the meter is in the concentrat ion mode.

3.

Adjust the meter to the concentration of the standard and fix the value in the memory according to the meter manuf acturer's instructions.

4. Rinse the electrodes with distilled water and blot dry.

To a 150 ml beaker, add 100 ml of the higher value s tandard and 1 ml of ISA. Place the beaker on a magnetic stirr er and begin stirring at a constant rate. Lower the electrode tips into the solution.

6.

Adjust the meter to the concentration of the standard and fix the value in the memory.

7.

Read the electrode slope according to the meter m anufacturer's instructions. Correct electrode operation is ind icated by a slope of 93-100%. See the **<u>TROUBLESHOOTING</u>** sections if the slope is not within this range.

MEASUREMENT

Measuring Hints

All samples and standards should be at the same temperature for pre cise measurement. A difference of $1^{\circ}C$ in temperature will result in about a 2% measurement error.

Constant, but not violent, stirring is necessary for accurate measu rement. Magnetic stirrers can generate sufficient heat to change t he solution temperature. To counteract this effect, place a piece of insulating material, such as styrofoam sheet, between the stirre r and the beaker.

Always rinse the electrodes with distilled water and blot dry betwe enmeasurements. Use a clean, dry tissue to prevent cross contamination.

When making low level measurements (below $1X10^{-5}M$), use only plastic lab-ware and cover the beaker with Parafilm to avoid loss of cyani de. When making high cyanide measurements (above $1X10^{-3}M$), samples should be diluted before measurements.

Use fresh standards for calibration.

Use 1 ml of ISA for each 100 ml of sample or standard.

Always check to see that the membrane is free from air bubbles afte r immersion into the standard or sample.

Sample Requirements

5.

All samples must be aqueous and not contain organics which can diss olve the epoxy electrode body and/or the cement bonding the sensing crystal to the electrode body. Infrequent measurements in solutio ns containing methanol, benzene, or acetonitrile are permitted. Hi ghly polar solvents slowly attack the epoxy body electrode. Please check with pHoenix Electrode Company before using the electrode in other solvents.

The temperature of the sample solutions and of the standard solutio ns should be the same and below $80^{\circ}C$.

Interferences should be absent. If they are present, use the proce dure found in the **Interferences** and **Electrode Response** sections to remove them.

Measurements above $1\times10^{-3}M$ should be done infrequently, as cyanide i on slowly erodes the membrane. (See section on **Electrode Life**.) It may be necessary to polish the membrane occasionally with polishin g paper or jeweller's rouge as the electrode is used. Samples should be diluted below $1\times10^{-3}M$ if possible.

Proper pH is ensured if ISA is used. The pH should be above 10 so that cyanide is present as CN^{-1} rather than as HCN in all standards and samples.

Units of Measurement

Cyanide ions can be measured in units of ppm, moles per liter, or a ny other convenient concentration unit. Table 1 indicates some con centration units and conversion factors.

TABLE 1: Concentration Unit Conversion Factors

ppm CN ⁻	moles/liter
26.00	1.0X10 ⁻³
10.00	3.8×10^{-4}
2.60	1.0×10^{-4}
1.00	3.8X10 ⁻⁵
0.26	1.0X10 ⁻⁵

MEASUREMENT PROCEDURE

Direct Measurement

Direct measurement is a simple procedure for measuring a large numb er of samples. A single meter reading is all that is required for each sample. The ionic strength of samples and standards should be made the same by adjustment with ISA. The temperature of both sam ple solutions and standard solutions should be the same.

Direct Measurement of Cyanide (using a pH/mV meter)

1.

By serial dilution, prepare 10^{-3} , 10^{-4} , and 10^{-5} M or 10, 1 and 0.1 ppm standards for the cyanide ion from the 0.01M ppm or 1000 ppm standards. Prepare standards with a composition s imilar to the samples if the samples have an ionic strength above 0.1M.

2.

Place 100 ml of the most dilute standard solution in a 1 50 ml plastic beaker. Place the beaker on the magnetic stirre r and begin stirring at a constant rate. Add 1 ml of 10M NaOH (ISA). After assuring that the meter is in the mV mode, lower the electrode tips into the solution. After the reading has stabilized, record the mV reading.

3.

Place 100 ml of the mid-range solution in a 150 ml plast ic beaker. Place the beaker on the magnetic stirrer and begin stirring at a constant rate. Add one ml of 10M NaOH (ISA). A fter rinsing the electrodes with distilled water, blot dry, an d lower the electrode tips into the solution. After the readi ng has stabilized, record the mV reading.

4.

Place 100 ml of the most concentrated standard solution in a 150 ml plastic beaker. Place the beaker on the magnetic stirrer and begin stirring at a constant rate. Add 1 ml of 10 M NaOH (ISA). After rinsing the electrodes with distilled wat er, blot dry, and lower the electrode tips into the solution. After the reading has stabilized, record the mV reading.

5.

Using the semi-logarithmic graph paper, plot the mV r eading (linear axis) against the concentration (log axis). Ex trapolate the curve down to about $1.0 \times 10^{-5} M$. A typical calibra tion curve can be found in Figure 1.

A calibration curve is constructed on semi-logarithmic paper when u sing a pH/mV meter in the millivolt mode. The measured electr ode potential in mV (linear axis) is plotted against the stand ard concentration (log axis). In the linear region of the curv e, only three standards are necessary to determine a calibrati on curve. In the non-linear region, additional points must be measured. The direct measurement procedures given are for the linear portion of the curve. The non-linear portion of the c urve requires the use of low level procedures. To a clean, dry, 150 ml plastic beaker, add 100 ml of th e sample and 1 ml of 10M NaOH (ISA). Place the beaker on the magnetic stirrer and begin stirring. Rinse the electrodes wit h distilled water, blot dry, and lower the electrode tips into the solution. When the reading has stabilized, record the mV reading. Using the calibration curve, determine the sample c oncentration.

7.

6.

The calibration should be checked every 2 hours. Assumin g no change in ambient temperature, place the electrode tips i n the mid-range standard. After the reading has stabilized, c ompare it to the original reading recorded in Step 3 above. A reading differing by more than 0.5 mV or a change in the ambi ent temperature will necessitate the repetition of Steps 2-5 a bove. A new calibration curve should be prepared daily.

Direct Measurement of Cyanide (using an ion meter)

1.

By serial dilution of the $1.0X10^{-2}M$ or 1000 ppm cyanide s tandard, prepare two standards whose concentration is near the expected sample concentration. Measure 100 ml of each standa rdintoindividual150 ml beakers and add1ml of 10MNaOH(ISA) to each.

2.

Place the more dilute solution on the magnetic stirrer and begin stirring at a constant rate. Assure that the meter is in the concentration mode.

3. Lower the electrode tips into the solution.

4.

Adjust the meter to the concentration of the cyanide sta ndard and fix the value in the memory according to the meter m anufacturer's instructions after stabilization of the reading.

5. Rinse the electrodes with distilled water and blot dry.

6.

Place the more concentrated solution on the magnetic sti rrer and begin stirring at a constant rate.

7. Lower the electrode tips into the solution.

8.

Adjust the meter to the concentration of the cyanide sta ndard and fix the value in the memory according to the meter m anufacturer's instructions after stabilization of the reading. For low level measurements (below $8 \times 10^{-6} M$ or 0.2 ppm), pl ace the rinsed, dried electrodes into a solution containing 10 0 ml of distilled water and 1 ml ISA. After stabilization, fi x the blank value in the meter according to the meter manufact urer's instruction.

10.

9.

After rinsing the electrodes and blotting dry, place the electrode tips into 100 ml of the sample and 1 ml of ISA. Af ter stabilization, read the concentration directly from the me ter display.

11.

The calibration should be checked every 2 hours. Assumin g no change in ambient temperature, place the electrode tips i n the first cyanide standard. After the reading has stabilize d, compare it to the original reading in Step 4 above. A read ing differing by more than 0.5 mV or a change in the ambient t emperature will necessitate the repetition of Step 2-10 above. The meter should be re-calibrated daily.

Low-Level Cyanide Measurements (using a pH/mV meter)

Use the following low level cyanide measurement procedure in the no n-linear portion of the calibration curve (below $8X10^{-6}M$ or 0.2 ppm). (See Figure 1.) A more accurate electrode indicator technique, s uch as titration, using a silver/sulfide ion electrode, may be pref erred below these levels.

1.

By serial dilution, prepare 100 ml of $1.0 \times 10^{-3} M$ or 10 p pm cyanide standard from the $1.0 \times 10^{-2} M$ or 1000 ppm cyanide standard.

2.

Using a 150 ml plastic beaker, add 100 ml of distilled w ater and 1 ml ISA. Place the beaker on the magnetic stirrer and begin stirring at a constant rate. Lower the electrode t ips into the solution. Assure that the meter is in the mV mode.

3.

Increments of the standard should be added to the beaker according to the steps outlined in Table 2 below. After the reading stabilizes, record the mV reading for each addition.

8

		Added	Concent	tratio	n
Step 1	Pipet	Volume (ml)	ppm	1	<u>N</u>
1	A	0.1		0.01	1.0X10 ⁻⁶
2	A	0.1		0.02	2.0X10 ⁻⁶
3	A	0.2		0.04	4.0X10 ⁻⁶
4	A	0.2		0.06	6.0X10 ⁻⁶
5	A	0.4		0.10	9.9X10 ⁻⁶
б	В	2.0		0.29	2.9X10 ⁻⁵
7	В	2.0		0.48	4.8X10 ⁻⁵

TABLE 2: Low Level Measurement Calibration Curve

Pipet A = 1 ml graduated pipet

Pipet B = 2 ml pipet

Solutions: additions of 1.0X10⁻³M or 10 ppm standard to 100 ml of distilled water and 1 ml ISA

4.

On semi-logarithmic paper, plot the concentration (log a xis) against the millivolt reading (linear axis) as in Figure 1. Reserve the final solution for checking the electrode each hour.

5.

To a 150 ml plastic beaker, add 100 ml of sample and 1 m l of ISA. Place the beaker on the magnetic stirrer and begin stirring at a constant rate. After rinsing the electrodes, bl ot dry and lower the electrode tips into the solution. After stabilization of the reading, read the mV potential and determ ine the concentration from the calibration curve. A new low l evel calibration curve should be prepared daily using fresh st andards.

ELECTRODE CHARACTERISTICS

Reproducibility

Electrode measurements reproducible to +2% can be obtained if the e lectrode is calibrated every hour. Factors such as temperature flu ctuations, drift, and noise limit reproducibility. Reproducibility is independent of concentration within the electrode's operating range.

Interferences

A layer of silver metal may form on the electrode surface in the pr esence of strongly reducing solutions, such as photographic develop er. Ions forming very insoluble silver salts will cause the electr odes to malfunction if present in solution at sufficiently high lev els to form a layer of silver salt on the membrane surface. Electr ode performance can be restored by polishing if the surface of the sensing element becomes contaminated. See the section **Electrode Re sponse** for proper polishing procedure. Solutions containing oxidizing agents such as Fe^{+3} , Cu^{+2} , and MnO_4^- will not affect electrode performance. All samples must be fr ee of mercury.

The maximum allowable ratio of interfering ions to cyanide ions is given in Table 3. The ratio is expressed as the ratio of the inter fering ion concentration in moles per liter to the sample cyanide c oncentration in moles per liter. Readings will be in error if this ratio is exceeded. Neither the accuracy of the measurement nor th e surface of the electrode membrane will be affected if the ratio i s less than that listed in the table.

TABLE 3: Maximum Allowable Ratio of Interfering Ion to Cyanide Ion

Interferences	Maximum Ratio
Cl-1	1X10°
Br^{-1}	5X10 ³
I^{-1}	1×10^{-1}
S^{-2}	must be absent

When using the cyanide ion electrode, an example of the use of Tabl e 3 follows:

What is the maximum level of bromide allowable in a sample whose cy anide concentration is $1X10^{-5}M$? Using Table 3, the maximum ratio is: $[Br^{-1}] = 5X10^{3}$ $[CN^{-1}]$

 $[Br^{-1}] = 5X10^{3} [CN^{-1}]$ = (5X10³)(1X10⁻⁵) [Br⁻¹] = 5X10⁻²M = maximum bromide concentration for no interference

Complexation

Total concentration, Ct, consists of free ions, Cf, and complexed o r bound ions, Cc, in solutions:

Ct = Cf + Cc

Since the electrode only responds to free ions, any complexing agen t in the solution reduces the measured concentration of ions.

Hydrogen ions and many metal ions form complexes with cyanide ions. The presence of any complexing agent lowers the measured concentr ation. Since the electrode measures only free cyanide ions, use of 10MNaOH(ISA) is essential, since it eliminates complexation by hydrogen.

EDTA can be used to break up cyanide complexes formed with many met al ions, including cadmium, copper, nickel, and zinc. To a sample w hose cyanide concentration is not more than 10 ppm, or about 1×10^{-3} M , add acetic acid to make the sample solution's pH~4. Add EDTA (disodium) to about 0.02M (or about 0.76 grams Na₂ EDTA per 100 ml sample.) <u>In a hood</u>, heat the mixture to about 50°C for about five minutes to speed up the decomplexation. After cooling the solution , add 10M NaOH (ISA) until the pH~13. The cyanide remains free lon g enough for concentration measurements to be made, since EDTA comp lexes of the metals break up very slowly.

This method will not work for silver, mercury, gold, or cobalt, sin ce they will bind the cyanide too strongly.

Temperature Influences

Samples and standards should be within $\pm 1^{\circ}$ C of each other, since el ectrode potentials are influenced by changes in temperature. A 1°C difference in temperature results in a 2% error at the 1.0×10^{-3} M con centration level. Because of the solubility equilibria on which th e electrode depends, the absolute potential of the reference electr ode (or reference portion of the combination electrode) changes slo wly with temperature. The slope of the electrode, as indicated by the factor "S" in the Nernst equation, also varies with temperature . Table 4 gives values for the "S" factor in the Nernst equation for the cyanide ion.

TABLE 4: Temperature vs. Values for the Electrode Slope

Temp.(°C)	"S"
0	54.20
10	56.18
20	58.17
25	59.16
30	60.15
40	62.13
50	64.11

If changes in temperature occur, the electrodes should be re-calibrated.

The temperature range for the pHoenix Cyanide Ion Electrodes is 0° - 80° C, provided that temperature equilibrium has occurred. If the temperature varies substantially from room temperature, equilibrium t imes up to one hour are recommended.

Electrode Response

Plotting the electrode mV potential against the cyanide concentrati on on semi-logarithmic paper results in a straight line with a slop e of about 57 mV per decade. (Refer to Figure 1.)

The time needed to reach 99% of the stable electrode potential read ing, the electrode response time, varies from several seconds in hi ghly concentrated solutions to several minutes near the detection l imit. (Refer to Figure 2.)

A drifting potential reading or a decrease in electrode slope may m ean that the electrode membrane needs polishing.

To polish the membrane:

1. If using polishing paper, cut off a 1-2" piece and place it face up on the lab bench.

2.

Put a few drops of distilled or deionized water in the c enter of the paper.

3.

Holding the paper (cotton) steady with one hand, bring t he membrane of the electrode down perpendicular to the paper and, with a slight swirling motion, gently polish the tip of t he electrode against the surface of the polishing paper (cotto n) for a few seconds.

4.

Rinse the electrode surface with distilled or deionized water and soak the electrode tip in standard solution for about t five minutes before use.

5.

If using jeweller's rouge, place a cotton ball on the ta ble top and flatten it using the bottom of a beaker.

б.

Put 1-2 drops of distilled or deionized water in the center of the cotton pad.

7. Add a small amount of jeweller's rouge to the damp cotton.

8. Continue with Steps 3 and 4 above.

Rinse the electrode surface with distilled water and soak the elect rode tip in $1X10^{-4}M$ or 1 ppm standard solution for about two minutes before use.

Limits of Detection

Cyanide levels from $5X10^{-6}M$ to $1X10^{-2}M$ cyanide can be measured with the cyanide electrodes. However, since cyanide ion attacks the ele ctrode membrane, measurements above $1X10^{-3}M$ should be done only inte rmittently.

The electrodes respond to cyanide in the sample as well as to ions dissolved from the membrane at low levels. The electrode membrane shows a very slight water solubility. The detection limit of the e lectrode is determined by this factor. Figure 1 shows the theoreti cal linear response -- the dashed line -- in comparison with the ac tual response (solid line). The low-level procedure is recommended if measurements are to be made in the non-linear region below $8 \times 10^{-6} M$.

Plastic lab-ware must be used and the beakers must be covered with Parafilm for low level cyanide determinations or cyanide will be lo st. Allow a longer stabilization time before taking the meter read ing for best results.

pH Effects

The cyanide electrode can be used over the pH range 11 to 13. It is necessary to adjust the sample pH to above 11 using the recommend ed ISA to convert all cyanic acid species in solution to cyanide.

Electrode Life

The cyanide electrode will last six months in normal laboratory use . On-line measurements might shorten operational lifetime to sever al months. In time, the response time will increase and the calibr ation slope will decrease to the point calibration is difficult and electrode replacement is required.

Electrode Storage

The cyanide electrode may be stored for short periods of time in 1. 0×10^{-4} M cyanide solution with ISA added. For longer storage (longer than two weeks), rinse and dry the sensing pellet and cover the me mbrane tip with any protective cap shipped with the electrode. The reference portion of the combination electrode (or the outer chamb er of the reference electrode) should be drained of filling solutio n, if refillable, and the rubber insert placed over the filling hole.

ELECTRODE THEORY

Electrode Operation

The pHoenix Cyanide Ion Electrodes consist of a solid membrane cont aining a mixture of inorganic silver compounds bonded into the tip of a glass or epoxy electrode body. An electrode potential develop s across the membrane when the electrode is in contact with solutio n containing cyanide ions and is capable of measuring free cyanide ions. This potential is measured against a constant reference pote ntial, using a standard pH/mV meter or an ion meter. The Nernst eq uation describes the level of cyanide ions in solution correspondin g to the measured potential:

 $E = E_{o} - S \log X$

where:

E = measured electrode potential

 E_{\circ} = reference potential (a constant)

S = electrode slope (~57 mV/decade)

X = level of cyanide ions in solution

The activity, X, represents the effective concentration of free cya nide ions in the solution. Total cyanide concentration, Ct, may in clude some bound as well as free cyanide ions. Since the electrode only responds to free ions, the concentration of the free ions, Cf , is found by:

Cf = Ct - Cb

where Cb represents the concentration of all bound or complexed cya nide ions.

The activity is related to the free ion concentration, C_f , by the a ctivity coefficient, a, by:

 $X = a C_f$

Activity coefficients vary, depending on total ionic strength, I, d efined as:

 $I = 1/2 O C_x Z_x^2$

where: C_x = concentration of ion X Z_x = charge of ion X O = sum of all of the types of ions in the solution

In the case of high and constant ionic strength relative to the sen sed ion concentration, the activity coefficient, a, is constant and the activity, X, is directly proportional to the concentration. All samples and standards containing cyanide ions have ionic streng th adjuster (ISA) added so that the background ionic strength is hi gh and constant relative to variable concentrations of cyanide. Th e recommended ISA for the cyanide electrode is sodium hydroxide, Na OH, though other basic solutions can be used as long as they do not contain ions that would interfere with the electrode's response to cyanide.

The reference electrode must also be considered. When two solution s of different composition are brought into contact with one anothe r, liquid junction potentials arise. Millivolt potentials occur fr om the inter-diffusion of ions in the two solutions. Electrode cha rge will be carried unequally across the solution boundary resultin g in a potential difference between the two solutions, since ions d iffuse at different rates. When making measurements, it is importa nt to remember that this potential be the same when the reference i s in the standardizing solution as well as in the sample solution o r the change in liquid junction potential will appear as an error i n the measured electrode potential.

The composition of the liquid junction filling solution in the refe rence electrode is most important. The speed with which the positi ve and negative ions in the filling solution diffuse into the sampl e should be equitransferent. No junction potential can result if t he rate at which positive and negative charge carried into the samp le is equal.

TROUBLESHOOTING GUIDE

The goal of troubleshooting is the isolation of a problem through c hecking each of the system components in turn: the meter, the plast ic-ware, the electrodes, the standards and reagents, the sample, an d the technique.

Meter

The meter may be checked by following the check-out procedure in th e instrument instruction manual.

Plastic-ware

Clean plastic-ware is essential for good measurement. Be sure to w ash the plastic-ware well with a mild detergent and rinse very well with distilled or deionized water.

Electrodes

The electrodes may be checked by using the procedure found in the s ections entitled **Electrode Slope Check**.

1.

Be sure to use distilled or deionized water when following the procedures given in **Electrode Slope Check**.

2.

If the electrode fails to respond as expected, see the sections **Measuring Hints** and **Electrode Response**. Repeat the s lope check.

3.

If the electrodes still fail to respond as expected, sub stitute another cyanide ion electrode that is known to be in g ood working order for the questionable electrode. If the prob lem persists and you are using an electrode pair, try the same routine with a working reference electrode.

4.

If the problem persists, the reagent may be of poor qual ity, interferences in the sample may be present or the techniq ue may be faulty. See **Reagents**, **Sample**, and **Technique** section s below.

5.

If another electrode is not available for test purposes, or if the electrode in use is suspect, review the instruction manual and be sure to:

- Clean and rinse the electrodes thoroughly.

- Prepare the electrodes properly.

- Use the proper filling solution.

- Adjust the pH and the ionic strength of the solution by t he use of the proper ISA.

- Measure correctly and accurately.

- Review **TROUBLESHOOTING** HINTS.

Standards & Reagents

Whenever problems arise with the measuring procedure that has been used successfully in the past, be sure to check the standard and re agent solutions. If in doubt about the credibility of any of the s olutions, prepare them again. Errors may result from contamination of the ISA, incorrect dilution of standards, poor quality distille d/deionized water, or a simple mathematical miscalculation.

Sample

Look for possible interferences, complexing agents, or substances w hich could affect the response or physically damage the sensing ele ctrode (or the reference electrode) if the electrodes work perfectl y in the standard, but not in the sample.

Try to determine the composition of the samples prior to testing to eliminate a problem before it starts. (See Measuring Hints, Sampl e Requirements, and Interferences.)

Technique

Be sure that the electrode's limit of detection has not been exceed ed. Be sure that the analysis method is clearly understood and is compatible with the sample.

Refer to the instruction manual again. Reread **<u>GENERAL</u> <u>PREPARATION</u>** and **ELECTRODE CHARACTERISTICS**.

If trouble still persists, call pHoenix Electrode Company at 1 -800-522-7920 and ask for the Technical Services Department.

TROUBLESHOOTING HINTS

Symptom	Possible Causes	Next Step
Out of Range Reading	defective meter (i	check meter with shorting strap see meter nstruction
	manual)	
	defective electrode	check electrode operation
plugged	electrodes not in properly and :	unplug electrodes reseat
	reference electrode not filled	be sure reference electrode is filled
	air bubble on membrane	remove bubble by re-dipping electrode
	electrodes not in solution	put electrodes in solution

Noisy or Unstable defective meter check meter with Readings (readings shorting strap continuously or air bubble on remove bubble by rapidly changing) membrane re-dipping electrode electrode exposed soak electrode in to interferences cyanide standard ISA not used use recommended ISA meter or stirrer ground meter or not grounded stirrer defective electrode replace electrode Drift (reading samples and standards allow solutions to slowly changing at different come to room temperin one direction) temperatures ature before measurement electrode exposed check section to complexing entitled Complexation agents incorrect reference use recommended filling solution filling solution Low Slope or standards contaminprepare fresh No Slope ated or incorrectly standards made ISA not used use recommended ISA standard used as ISA use ISA electrode exposed check section entitled to complexing agents Complexation air bubble on remove bubble by membrane re-dipping probe "Incorrect Answer" plot millivolts on incorrect scaling (but calibration of semi-log paper the linear axis. On curve is good) the log axis, be sure concentration numbers within each decade are

	increasing with increasing concentration.
incorrect sign	be sure to note sign of millivolt number correctly
incorrect standard	ls prepare fresh standards
wrong units used	apply correct conversion factor: 10 ⁻³ M = 26.0 ppm CN ⁻
complexing agents in sample	check section entitled Complexation

SPECIFICATIONS

Concentration Range:	1X10 ⁻² M to 5X10 ⁻⁶ M 260 to 0.1 ppm
pH Range:	11-13
Temperature Range:	0°-80°C
Resistance:	<1 Mohms
Reproducibility:	<u>+</u> 2%
Size:	110 mm length 12 mm diameter 1 m cable length
Storage:	store in cyanide standard with ISA added

ORDERING INFORMATION

P/N	DESCRIPTION
CN01501	Cyanide Ion Electrode, mono (reference electrode necessary), epoxy body
CN01502	Cyanide Ion Electrode, combination, glass body
CN01503	Cyanide Ion Electrode, combination, epoxy body
5731429	Reference Electrode, double junction, sleeve, epoxy body, for use with the CN01501
CN0IS01	Ionic Strength Adjuster (ISA), 10M NaOH
R001011	Electrode Filling Solution, 4M KCl/AgCl, 125 ml, for the 5731429 inner chamber
R001015	Electrode Filling Solution, 1M $\rm KNO_3$, 125 ml, for the 5731429 outer chamber and the CN01502 electrode
R001046	Electrode Filling Solution, 1M $\rm KNO_3/0.1M$ KCl, 125 ml, for the CN01503 electrode
948201	Polishing Paper for the Cyanide Ion Electrodes

TABLE OF CONTENTS

General Instructions
General Preparation
Measurement
Measurement Procedure
<pre>(using an ion meter)7 low-level cyanide measurements (using a pH/mV meter)</pre>
Electrode Characteristics9reproducibility9interferences9complexation10temperature influences11electrode response12limits of detection13pH effects13electrode life13electrode storage14
Electrode Theory
Troubleshooting Guide15
Troubleshooting Hints17
Specifications
Ordering Information
Table of Contents