

# Calibration Station

Specification Sheet (Part/REF # 243575)



Calibration station used for validation / calibration of sensors in gaseous media. Equipped with two inlets for a wide array of test scenarios.

## Product Specifications

<b>For Sensors with</b>	120 mm, 12 mm, PG13,5
<b>Dimensions (W x D x H)</b>	137 x 115 x 186 mm
<b>Housing Material</b>	Black PEEK with aluminium holder
<b>Ø O.D. Fitting</b>	8,7 mm
<b>Ø I.D. Fitting</b>	3,5 mm
<b>Pressure Range bar g</b>	0 to 12 bar

Specifications are subject to change without notice

**Spec. Version B**

<https://www.hamiltoncompany.com/process-analytics/sensors/243575>